



# Waters

THE SCIENCE OF WHAT'S POSSIBLE.™



The Pre-Slit Caps are manufacturer by making a slit in the silicone layer, but leaving the PTFE barrier uncut

**IT TECH RESEARCH (M) SDN BHD**  
No 22(A), Ground Floor,  
Jalan Putra Mahkota 7/8A,  
Section 7, Putra Heights,  
47650 Subang Jaya, Selangor  
Tel: 03-5192 7742 Fax: 03-5192 7743  
Email: sales@ittech.com.my

# PROMOTION

Purchase **Waters Bonded Pre-Slit Cap (RM180)**,  
To get Free Gift of your choice  
(contact us for more details)



**PROMO CODE : IT03**

Validity : While Stock Last

\*Terms & Conditions Applied



- Less pressure required for autosampler to penetrate septa
- Avoiding bent or broken needles
- Provide a larger target area for needle to penetrate even when needle doesn't come down in center of septa
- Uncut PTFE layer provide better seal compared to other slitted septa