

PETTE vario

/// Simple liquid handling for variable volumes

An ergonomic feel combined with uncompromising precise performance: the IKA vario pipettes can be used for many applications in the modern laboratory. Colour coding simplifies the rapid selection of a suitable pipette. Resistance to impact, UV and chemicals ensures that the tip cone provides long-term reliable dosing. The PETTE vario are fully autoclavable, therefore ensuring easy sterilization. This guarantees fault-free, precise and safe operation.



Features

Ergonomic design

- Exchangeable grips in different shapes and materials with minimal frictional and spring forces
- One-handed adjustment of the volume using the multifunctional knob.



Durable construction

The tip cone and piston are coated with a high-quality diamond-like carbon (DLC) that guarantees particularly robust pipettes and reliably meter in the long term.



Special accs (5ml/10ml only)

- Filter: to protect pipette/sample from contamination during sample suck up
- Grease: as lubricant



A large display

- the set volumes in easily reads characters.

User-friendly calibration

- Software used for calibration tests that are user-friendly.

Simple maintenance and adjustment

- The volume unit with piston and seal can be disassembled and cleaned in few steps.
- If adjustment/calibration is necessary, it can be done without any additional tool.

***Universal tips usage for all brand of tips.**

Overview PETTE vario

PETTE VARIO	P/N	MAX. VOLUME TIP
0.1 – 2 µl	0020011210	10 µl
0.5 – 10 µl	0020011211**	10 µl
2 – 20 µl	0020011213**	200 µl
10 – 100 µl	0020011214**	200 µl
20 – 200 µl	0020011215**	200 µl
100 – 1,000 µl	0020011216**	1,000 µl
0.5 – 5 ml	0020011217	5 ml
1 – 10 ml	0020011218	10 ml

PRICE NOT MORE THAN RM600



EXSTOCK AVAILABLE.

** Free micropipette tips with purchase of each micropipette