

GAS PURIFIERS AND CONSUMABLES

Scientific Glass Technology is the manufacturer of a wide range of gas purifiers and related consumables used in the chromatography market. It also manufactures and sell various kinds of laboratory tools such as capillary column cutters and many more.

GAS PURIFICATION SYSTEM



Super Clean Gas Purification systems make changing gas purifiers quick and easy. A base plate allows purifiers to be exchanged without introducing oxygen.

INLINE GAS PURIFICATION SYSTEM



Click-On Inline Super-Clean Big Traps are inline traps designed with Click-On adaptor connectors which allows inline cartridges to be exchanged without introducing oxygen.

ELECTRONIC MAINTENANCE INDICATOR



The Electronic Indicator Device progressively displays real-time consumable maintenance status instead of a 'last-minute' colour change from traditional indicators.

IN-LINE BIG TRAP GAS PURIFICATION SYSTEM



Click-On Inline Super-Clean Big Traps are inline traps designed with Click-On adaptor connectors which allows inline cartridges to be exchanged without introducing oxygen.

LABORATORY TOOLS



The Capillary Column Cutter for GC is especially designed to allow everybody, independent from skill and experience level, to make the perfect cut.

Quality
Credentials



Management
System
ISO 9001:2015

www.tuv.com
ID 9108636629

Get back to us for the best pricing.

For more info, you may visit
<https://www.sgt-nl.com/Default.aspx>